



Uni - Layer Shade Selection

Shade	W0	W1	W2	W3	A1	A2	A3	A3.5	A4	Root-A	B1	B2	B3	B4	C1	C2	C3	C4	D2	D3	D4	
1st Opaque	WoO	W1O	W2O	W3O	A1O	A2O	Base	Base	Base	Base	B1O	B2O	Base	Base	C1O	C2O	Base	Base	Base	Base	Base	
2nd Opaque	W0O	W1O	W2O	W3O	A1O	A2O	A3O	A3.5	A4O	root AO	B1O	B2O	B3O	B4O	C1O	C2O	C3O	C4O	D2O	D3O	D4O	
Crowns Anterior																						
1 – 2 mm	LT	LT	LT	LT	LT	LT	LT	LT	HT	HT	HT	LT	LT	LT	HT	LT	LT	LT	HT	LT	LT	HT
Crowns Posterior																						
1 – 2 mm	LT	LT	LT	LT	LT	LT	LT	LT	HT	HT	HT	LT	LT	LT	HT	LT	LT	LT	HT	LT	LT	HT
Bridges																						
1 – 2 mm	LT	LT	LT	LT	LT	LT	LT	LT	HT	HT	HT	LT	LT	LT	HT	LT	LT	LT	HT	LT	LT	HT
Pontic	OD-WO	OD-WO:2 OD-N :1	OD-WO:1 OD-N :2	OD-NM	OD-A1	OD-A2	OD-A3	OD-A3.5	OD-A4	OD-root A	OD-N :1 OD-B2 :1	OD-B2	OD-B2:1 OD-B4:1	OD-B4	OD-N :1 OD-C2 :1	OD-C2	OD-C2:1 OD-C4:1	OD-C4	OD-N :1 OD-D3 :1	OD-D3	OD-D3:1 OD-B4:1	
Art Stain																						
Cervical					A-Shade	A-Shade	A-Shade	A-Shade	A-Shade	A-Shade	B-Shade	B-Shade	B-Shade	B-Shade	C-Shade	C-Shade	C-Shade	C-Shade	D-Shade	D-Shade	D-Shade	
Incisal / Occlusal	Bl.-Gray White O-Brown	Bl.-Gray White O-Brown	Bl.-Gray White O-Brown	Bl.-Gray White O-Brown	Bl.-Gray White O-Brown	Bl.-Gray White O-Brown	Bl.-Gray White O-Brown	Bl.-Gray White O-Brown	Bl.-Gray White O-Brown	Bl.-Gray White O-Brown	Bl.-Gray White O-Brown	Bl.-Gray White O-Brown	Bl.-Gray White O-Brown	Bl.-Gray White O-Brown	Bl.-Gray White O-Brown	Bl.-Gray White O-Brown	Bl.-Gray White O-Brown	Bl.-Gray White O-Brown	Bl.-Gray White O-Brown	Bl.-Gray White O-Brown	Bl.-Gray White O-Brown	

SHOFU DENTAL GMBH
 Am Brüll 17
 40878 Ratingen
www.shofu.de