



VINTAGE AL & ZR

COLOR TABLE

Shade	W1	W2	W3	A1	A2	A3	A3.5	A4	rootA
Opaque Liner									
	W1O	W2O	W3O	A1O	A2O	A3O	A3.5O	A4O	rootAO
Cervical	—	—	—	—	—	A3B:2 AC:1	A3.5B:1 AC:1	A4B:1 AC:1	
	—	—	—	—	—	A3B:2 AC:1	A3.5B:1 AC:1	A4B:1 AC:1	AC
Margin		W1M:2 NM:1	W1M:1 NM:2		NM:1 A3M:1		A3M:2 rootAM:1	A3M:1 rootAM:2	
	W1M	W1M:2 NM:1	W1M:1 NM:2	NM	NM:1 A3M:1	A3M	A3M:2 rootAM:1	A3M:1 rootAM:2	rootAM
Opaque Dentin		OD-W1:2 OD-N:1	OD-W1:1 OD-N:2		OD-N:1 OD-A3:1		OD-A3:2 OD-rootA:1	OD-A3:1 OD-rootA:2	
	OD-W1	OD-W1:2 OD-N:1	OD-W1:1 OD-N:2	OD-N	OD-N:1 OD-A3:1	OD-A3	OD-A3:2 OD-rootA:1	OD-A3:1 OD-rootA:2	OD-rootA
Body									
	W1B	W2B	W3B	A1B	A2B	A3B	A3.5B	A4B	rootAB
Enamel		Opal 56:2 Opal 57:1	Opal 56:1 Opal 57:2				Opal 59:1 Opal 60:1		—
	Opal 56	Opal 56:2 Opal 57:1	Opal 56:1 Opal 57:2	Opal 57	Opal 58	Opal 59	Opal 59:1 Opal 60:1	Opal 60	—

Shade	B1	B2	B3	B4	C1	C2	C3	C4	D2	D3	D4
Opaque Liner											
	B1O	B2O	B3O	B4O	C1O	C2O	C3O	C4O	D2O	D3O	D4O
Cervical	—	—	B3B:1 BC:1		—	C2B:2 CC:1	C3B:1 CC:1				
	—	—	B3B:1 BC:1	BC	—	C2B:2 CC:1	C3B:1 CC:1	CC	DC		
Margin	NM:1 B2M:1		B2M:1 B4M:1		NM:1 C2M:1		C2M:1 C4M:1		NM:1 D3M:1		D3M:1 B4M:1
	NM:1 B2M:1	B2M	B2M:1 B4M:1	B4M	NM:1 C2M:1	C2M	C2M:1 C4M:1	C4M	NM:1 D3M:1	D3M	D3M:1 B4M:1
Opaque Dentin	OD-N:1 OD-B2:1		OD-B2:1 OD-B4:1		OD-N:1 OD-C2:1		OD-C2:1 OD-C4:1		OD-N:1 OD-D3:1		OD-D3:1 OD-B4:1
	OD-N:1 OD-B2:1	OD-B2	OD-B2:1 OD-B4:1	OD-B4	OD-N:1 OD-C2:1	OD-C2	OD-C2:1 OD-C4:1	OD-C4	OD-N:1 OD-D3:1	OD-D3	OD-D3:1 OD-B4:1
Body											
	B1B	B2B	B3B	B4B	C1B	C2B	C3B	C4B	D2B	D3B	D4B
Enamel										—	
	Opal 57	Opal 58	Opal 59	Opal 60	Opal 58	Opal 59	Opal 60	Opal 58	Opal 59	—	Opal 59

Opal Effect								
	Opal T	Opal SL	Opal WE	Opal MI	Opal OC	Opal AM-R	Opal AM-Y	Opal AM-V

Enamel Effect						Cervical Trans					
	BT	OT	PT	GT	T-GLASS	Cervical Trans	CT-CL	CT-W	CT-A	CT-B	CT-R

Color Effect									Gum Color		
	MP	MY	Mlv	RED	Y	O	G	W	Gum Color	Gum-LP	Gum-DP
						Correction					
	OM-Y	OM-LP	OM-DP	LPM	CLM	Correction			ADD-ON B	ADD-ON T	CPM-FINE



VINTAGE AL & ZR



Opal MI



Opal WE



Opal OC



Opal 56



Opal 57



Opal 58



Opal 59



Opal 60



Opal SL



Opal T



Opal AM-Y



Opal AM-V



Opal AM-R



OT



PT



BT



GT



T-GLASS



MY



Miv



MP



Y



O



RED



G



W



CT-CL



CT-W



CT-A



CT-B



CT-R